



TROVE2 Antibody, Biotin conjugated

Product Code	CSB-PA13749D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10155
Immunogen	Recombinant Human 60 kDa SS-A/Ro ribonucleoprotein protein (3-535AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60 kDa SS-A/Ro ribonucleoprotein (60 kDa Ro protein) (60 kDa ribonucleoprotein Ro) (RoRNP) (Ro 60 kDa autoantigen) (Sjogren syndrome antigen A2) (Sjogren syndrome type A antigen) (SS-A) (TROVE domain family member 2), TROVE2, RO60 SSA2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TROVE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.