




HLA-DQA1 Antibody, Biotin conjugated

Product Code	CSB-PA14839D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01909
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DQ alpha 1 chain protein (24-213AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DQ alpha 1 chain (DC-1 alpha chain) (DC-alpha) (HLA-DCA) (MHC class II DQA1), HLA-DQA1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DQA1

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA14839D0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.