




# eutC Antibody

<b>Product Code</b>	CSB-PA18239A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P19636
<b>Immunogen</b>	Recombinant Escherichia coli Ethanolamine ammonia-lyase light chain protein (1-295AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Escherichia coli
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ethanolamine ammonia-lyase light chain (EC 4.3.1.7) (Ethanolamine ammonia-lyase small subunit), eutC
<b>Immunogen Species</b>	Escherichia coli
<b>Research Area</b>	Others
<b>Target Names</b>	eutC

## Image



### Western Blot

Positive WB detected in: Rosseta bacteria lysate at 80µg, 40µg, 20µg, 10µg

All lanes: eutC antibody at 3µg/ml

### Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution  
predicted band size: 32 kDa

observed band size: 32 kDa