



glpK Antibody, FITC conjugated

Product Code	CSB-PA18407C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A6F3
Immunogen	Recombinant Escherichia coli Glycerol kinase protein (2-502AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycerol kinase (EC 2.7.1.30) (ATP:glycerol 3-phosphotransferase) (Glycerokinase) (GK), glpK
Immunogen Species	Escherichia coli (strain K12)
Research Area	Signal Transduction
Target Names	glpK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.