




Mmp13 Antibody

Product Code	CSB-PA17739A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33435
Immunogen	Recombinant Mouse Collagenase 3 protein (105-472AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagenase 3 (EC 3.4.24.-) (Matrix metalloproteinase-13) (MMP-13), Mmp13
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Mmp13

Image



IHC image of CSB-PA17739A0Rb diluted at 1:200 and staining in paraffin-embedded human colon cancer performed on a Leica Bond™ system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.