



Mmp3 Antibody, FITC conjugated

Product Code	CSB-PA17509C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28862
Immunogen	Recombinant Mouse Stromelysin-1 protein (104-477AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stromelysin-1 (SL-1) (EC 3.4.24.17) (EMS-2) (Matrix metalloproteinase-3) (MMP-3) (Transin-1), Mmp3
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Mmp3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.