



Ofloxacin Monoclonal Antibody

Product Code	CSB-MA003311I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Ofloxacin-BSA conjugate
Raised In	mouse
Specificity	no significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95% protein G purified
Isotype	IgG2b
Clonality	monoclonal
Alias	Oflocee, oflocet, oflocin, orocin, tarivid
Product Type	Small Molecules Antibodies
Target Names	OFX
Accession NO.	7C1F3