



# Sulfadimidine Monoclonal Antibody

<b>Product Code</b>	CSB-MA004011I0m
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Sulfadimidine-BSA conjugate
<b>Raised In</b>	mouse
<b>Specificity</b>	no significant cross-reactivity or interference was observed
<b>Tested Applications</b>	ELISA
<b>Form</b>	liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Caprylic Acid Ammonium Sulfate Precipitation purified
<b>Isotype</b>	IgG2a
<b>Clonality</b>	monoclonal
<b>Alias</b>	SM2, S-Mez, Sulfamethazine,4,6-Dimethylsulfadiazine, 4-Amino-N-(4,6-dimethyl-2-pyrimidinyl)benzenesulf-onamide
<b>Product Type</b>	Small Molecules Antibodies
<b>Target Names</b>	SM2
<b>Accession NO.</b>	2C1D2