



Sulfadimidine Monoclonal Antibody

Product Code	CSB-MA004011I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Sulfadimidine-BSA conjugate
Raised In	mouse
Specificity	no significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Caprylic Acid Ammonium Sulfate Precipitation purified
Isotype	IgG2a
Clonality	monoclonal
Alias	SM2, S-Mez, Sulfamethazine, 4,6-Dimethylsulfadiazine, 4-Amino-N-(4,6-dimethyl-2-pyrimidinyl)benzenesulf-onamide
Product Type	Small Molecules Antibodies
Target Names	SM2
Accession NO.	2C1D2