



Sulfametoxydiazine Monoclonal Antibody

Product Code	CSB-MA004111I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Sulfametoxydiazine-BSA conjugate
Raised In	mouse
Specificity	No significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%?protein G purified
Isotype	IgG1
Clonality	monoclonal
Alias	SMD, aristebon, sulfametoxydiazine, sulfameter,bayrena,durenat, 5-Methoxysulfadiazine
Product Type	Small Molecules Antibodies
Target Names	SMD
Accession NO.	5B11G3