



Chenodeoxycholic acid Monoclonal Antibody

Product Code	CSB-MA005711I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Chenodeoxycholic acid-BSA conjugate
Raised In	mouse
Specificity	no significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95% protein G purified
Isotype	IgG2b
Clonality	monoclonal
Alias	Chenodiol, chenix, chenic acid, chenodeoxycholate, gallodesoxycholic acid, chendol, chenofalk, chenodesoxycholic acid, CDCA
Product Type	Small Molecules Antibodies
Target Names	CDCA
Accession NO.	4E1D4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.