



Chenodeoxycholic acid Monoclonal Antibody

Product Code	CSB-MA005711I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Chenodeoxycholic acid-BSA conjugate
Raised In	mouse
Specificity	no significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%?protein G purified
Isotype	IgG2b
Clonality	monoclonal
Alias	Chenodiol, chenix, chenic acid, chenodeoxycholate, gallodesoxycholic acid, chendol, chenofalk, chenodesoxycholic acid, CDCA
Product Type	Small Molecules Antibodies
Target Names	CDCA
Accession NO.	4E1D4