



Hist1h2bm Antibody, FITC conjugated

Product Code	CSB-PA010413LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10854
Immunogen	Recombinant Mouse Histone H2B type 1-M protein (2-126AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H2B type 1-M (H2B 291B), Hist1h2bm
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Hist1h2bm
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.