



PTK2 Antibody

Product Code	CSB-PA14445A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00944
Immunogen	Recombinant Chicken Focal adhesion kinase 1 protein (376-683AA)
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Focal adhesion kinase 1 (FADK 1) (EC 2.7.10.2) (Focal adhesion kinase-related nonkinase) (FRNK) (p41/p43FRNK) (Protein-tyrosine kinase 2) (p125FAK) (pp125FAK), PTK2, FAK FAK1
Immunogen Species	Chicken
Research Area	Others
Target Names	PTK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.