



SYT1 Antibody

Product Code	CSB-PA023029LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21579
Immunogen	Recombinant Human Synaptotagmin-1 protein (99-416AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptotagmin-1 (Synaptotagmin I) (Sytl) (p65), SYT1, SVP65 SYT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SYT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.