



GAPDH Antibody, FITC conjugated

Product Code	CSB-PA009232HC01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00356
Immunogen	Recombinant Chicken Glyceraldehyde-3-phosphate dehydrogenase protein (2-333AA)
Raised In	Rabbit
Species Reactivity	Chicken
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glyceraldehyde-3-phosphate dehydrogenase (GAPDH) (EC 1.2.1.12) (Peptidyl-cysteine S-nitrosylase GAPDH) (EC 2.6.99.-), GAPDH, GAPD
Immunogen Species	Chicken
Research Area	Others
Target Names	GAPDH