



Bcl2 Antibody, FITC conjugated

Product Code	CSB-PA002611LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10417
Immunogen	Recombinant Mouse Apoptosis regulator Bcl-2 protein (5-205AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptosis regulator Bcl-2, Bcl2, Bcl-2
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Bcl2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.