




Bcl2 Antibody

Product Code	CSB-PA002611LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10417
Immunogen	Recombinant Mouse Apoptosis regulator Bcl-2 protein (5-205AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptosis regulator Bcl-2, Bcl2, Bcl-2
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Bcl2

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA002611LA01MO at dilution of 1:20

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.