




# HDAC4 Antibody

<b>Product Code</b>	CSB-PA17884A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P56524
<b>Immunogen</b>	Recombinant Human Histone deacetylase 4 protein (1-227AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Histone deacetylase 4 (HD4) (EC 3.5.1.98), HDAC4, KIAA0288
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	HDAC4

## Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA17884A0Rb at dilution of 1:20

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.