




STX6 Antibody

Product Code	CSB-PA022902LA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5ZL19
Immunogen	Recombinant Chicken Syntaxin-6 protein (4-231AA)
Raised In	Rabbit
Species Reactivity	Chicken, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-6, STX6
Immunogen Species	Chicken
Research Area	Others
Target Names	STX6

Image



Western blot

All lanes: Syntaxin-6 antibody at 2µg/ml

Lane 1: Mouse brain tissue

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 29 kDa

Observed band size: 29 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.