



PRDX1 Antibody, Biotin conjugated

Product Code	CSB-PA01165D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06830
Immunogen	Recombinant Human Peroxiredoxin-1 protein (1-199AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxiredoxin-1 (EC 1.11.1.15) (Natural killer cell-enhancing factor A) (NKEF-A) (Proliferation-associated gene protein) (PAG) (Thioredoxin peroxidase 2) (Thioredoxin-dependent peroxide reductase 2), PRDX1, PAGA PAGB TDPX2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRDX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.