



UBE2I Antibody, HRP conjugated

Product Code	CSB-PA00474B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63279
Immunogen	Recombinant Human SUMO-conjugating enzyme UBC9 protein (1-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SUMO-conjugating enzyme UBC9 (EC 2.3.2.-) (RING-type E3 SUMO transferase UBC9) (SUMO-protein ligase) (Ubiquitin carrier protein 9) (Ubiquitin carrier protein I) (Ubiquitin-conjugating enzyme E2 I) (Ubiquitin-protein ligase I) (p18), UBE2I, UBC9 UBCE9
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UBE2I
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.