



P4HB Antibody, FITC conjugated

Product Code	CSB-PA00254C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07237
Immunogen	Recombinant Human Protein disulfide-isomerase protein (29-508AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein disulfide-isomerase (PDI) (EC 5.3.4.1) (Cellular thyroid hormone-binding protein) (Prolyl 4-hydroxylase subunit beta) (p55), P4HB, ERBA2L PDI PDIA1 PO4DB
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	P4HB