




ARL2 Antibody

Product Code	CSB-PA01265A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36404
Immunogen	Recombinant Human ADP-ribosylation factor-like protein 2 protein (19-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor-like protein 2, ARL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARL2

Image



Western blot

All lanes: ARL2 antibody at 2µg/ml

Lane 1: EC109 whole cell lysate

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 21, 19 kDa

Observed band size: 21 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.