



MPZ Antibody, Biotin conjugated

Product Code	CSB-PA07844D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25189
Immunogen	Recombinant Human Myelin protein P0 protein (30-156AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin protein P0 (Myelin peripheral protein) (MPP) (Myelin protein zero), MPZ
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MPZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.