



MPZ Antibody, Biotin conjugated

Product Code	CSB-PA07844D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25189
Immunogen	Recombinant Human Myelin protein P0 protein (30-156AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin protein P0 (Myelin peripheral protein) (MPP) (Myelin protein zero), MPZ
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MPZ