



PSMB4 Antibody, FITC conjugated

Product Code	CSB-PA01374C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28070
Immunogen	Recombinant Human Proteasome subunit beta type-4 protein (46-264AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit beta type-4 (EC 3.4.25.1) (26 kDa prosomal protein) (HsBPROS26) (PROS-26) (Macropain beta chain) (Multicatalytic endopeptidase complex beta chain) (Proteasome beta chain) (Proteasome chain 3) (HsN3), PSMB4, PROS26
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMB4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.