




Tpt1 Antibody

Product Code	CSB-PA024134HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63028
Immunogen	Recombinant Mouse Translationally-controlled tumor protein (1-172AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Translationally-controlled tumor protein (TCTP) (21 kDa polypeptide) (p21) (p23), Tpt1, Trt
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Tpt1

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA024134HA01MO at dilution of 1:20

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.