



LRBA Antibody, Biotin conjugated

Product Code	CSB-PA013070LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50851
Immunogen	Recombinant Human Lipopolysaccharide-responsive and beige-like anchor protein (2-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipopolysaccharide-responsive and beige-like anchor protein (Beige-like protein) (CDC4-like protein), LRBA, BGL CDC4L LBA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LRBA