



Eno1 Antibody, Biotin conjugated

Product Code	CSB-PA007670HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17182
Immunogen	Recombinant Mouse Alpha-enolase protein (2-433AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-enolase (EC 4.2.1.11) (2-phospho-D-glycerate hydro-lyase) (Enolase 1) (Non-neural enolase) (NNE), Eno1, Eno-1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Eno1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.