




Eno1 Antibody, HRP conjugated

Product Code	CSB-PA007670HB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17182
Immunogen	Recombinant Mouse Alpha-enolase protein (2-433AA)
Raised In	Rabbit
Species Reactivity	Mouse, Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-enolase (EC 4.2.1.11) (2-phospho-D-glycerate hydro-lyase) (Enolase 1) (Non-neural enolase) (NNE), Eno1, Eno-1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Eno1

Image



Western Blot

Positive WB detected in: PC-3 whole cell lysate, HepG2 whole cell lysate, NIH/3T3 whole cell lysate, Rat heart tissue, Rat brain tissue, Rat skeletal muscle tissue, Mouse kidney tissue, Mouse liver tissue

All lanes: Eno1 antibody; HRP conjugated at 0.42µg/ml

Predicted band size: 48 kDa

Observed band size: 48, 34 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.