



AFP Antibody, Biotin conjugated

Product Code	CSB-PA09987D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02771
Immunogen	Recombinant Human Alpha-fetoprotein protein (19-609AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Binds copper, nickel, and fatty acids as well as, and bilirubin less well than, serum albumin. Only a small percentage (less than 2%) of the human AFP shows estrogen-binding properties.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-fetoprotein (Alpha-1-fetoprotein) (Alpha-fetoglobulin), AFP, HPAFP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	AFP