



Recombinant Bordetella petrii Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-YP018033BPB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A9IM91
Product Type	Recombinant Protein
Immunogen Species	Bordetella petrii (strain ATCC BAA-461 / DSM 12804 / CCUG 43448)
Purity	>85% (SDS-PAGE)
Sequence	MLFKDQLFLA SQYAAPHHLV SRLMGLAADC RTPAIKNWMI SRFVRRYGVD MSEALVQDPL AYDTFNQFFT RKL RADARPL DTEPGAVLCP ADGAISQLGP IEHGRIFQAK GHSYSLTSLG GGDPARAEPF MGGDFATIYL SPRDYHRVHM PCAGTLREMV HVPGR LFSVN PLTATHVPEL FARNERVA CL FDTEYGP MAM VLVGAMIVAS IETAWAGLVT PHKRQIRAQR YDEAARAPIH LDKGAEMGLF KLG
Source	Yeast
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.