




SPINK1 Antibody

Product Code	CSB-PA022575KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00995
Immunogen	Recombinant protein of Human SPINK1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;WB:1:500-1:2000
Storage Buffer	Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Affinity purification
Isotype	IgG
Alias	SPINK1;PCTT;PSTI;Spink3;TATI;TCP
Product Type	Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Intended Use	For research use only. Not for human, diagnostic or therapeutic use.
Target Names	SPINK1

Image



Western blot analysis of GST-SPINK1 recombinant protein using SPINK1 antibody.

Target Details

This protein is a trypsin inhibitor, which is secreted from pancreatic acinar cells into pancreatic juice. It is thought to function in the prevention of trypsin-catalyzed premature activation of zymogens within the pancreas and the pancreatic duct. Mutations in this gene are associated with hereditary pancreatitis and tropical calcific pancreatitis.