




SPINK1 Antibody

Product Code	CSB-PA022575KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00995
Immunogen	Recombinant protein of Human SPINK1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; WB: 1:500-1:2000
Form	Liquid
Storage Buffer	Store at -20oC or -80oC. Avoid freeze / thaw cycles. Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
Purification Method	Affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SPINK1

Image



Western blot analysis of GST-SPINK1 recombinant protein using SPINK1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.