



# Recombinant Human Nucleoside diphosphate-linked moiety X motif 6 (NUDT6)

<b>Product Code</b>	CSB-EP016176HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53370
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRQPLSWGRW RAMLARTYGP GPSAGYRWAS GAQGYVRNPP VGACDLQGEL DRFGGISVRL ARLDALDRLD AAFQKGLQA AVQQWRSEGR TAVWLHIPIL QSRFIAPAAS LGFCFHHAES DSSTLTLWLR EGPSRLPGYA SHQVGVAGAV FDESTRKILV VQDRNKLNK WKFPGLSEP EEDIGDTAVR EVFEETGIKS EFRSVLSIRQ QHTNPGAFGK SDMYIICRLK PYSFTINFCQ EECLRCEWMD LNDLAKTENT TPITSRVARL LLYGYREGFD KIDLTVEELP AVYTGLFYKL YHKELPENYK TMKGID
<b>Source</b>	E.coli
<b>Target Names</b>	NUDT6
<b>Protein Names</b>	Recommended name: Nucleoside diphosphate-linked moiety X motif 6 Short name= Nudix motif 6 Alternative name(s): Antisense basic fibroblast growth factor Protein GFG
<b>Expression Region</b>	1-316
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.