




COL5A1 Antibody

Product Code	CSB-PA005748KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20908
Immunogen	A synthetic peptide of Human COL5A1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Storage Buffer	Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Affinity purification
Isotype	IgG
Alias	COL5A1; collagen, type V, alpha 1; collagen alpha-1(V) chain
Immunogen Species	Homo sapiens (Human)
Intended Use	For research use only. Not for human, diagnostic or therapeutic use.
Target Names	COL5A1

Image



Western blot analysis of extracts of Lung cell line, using COL5A1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.