




CD83 Antibody

Product Code	CSB-PA004962KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01151
Immunogen	Recombinant protein of Human CD83
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Storage Buffer	Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Affinity purification
Isotype	IgG
Alias	CD83; BL11; HB15
Immunogen Species	Homo sapiens (Human)
Intended Use	For research use only. Not for human, diagnostic or therapeutic use.
Target Names	CD83

Image



Western blot analysis of extracts of K562 cell lines, using D83 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.