



# IGF2 Antibody

<b>Product Code</b>	CSB-PA011088KA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P01344
<b>Immunogen</b>	Recombinant protein of Human IGF2
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB; WB: 1:500-1:2000
<b>Form</b>	Liquid
<b>Storage Buffer</b>	Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
<b>Purification Method</b>	Affinity purification
<b>Isotype</b>	IgG
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Metabolism
<b>Target Names</b>	IGF2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.