




# SMARCA4 Antibody

<b>Product Code</b>	CSB-PA021801KA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P51532
<b>Immunogen</b>	A synthetic peptide of Human SMARCA4
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
<b>Form</b>	Liquid
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Affinity purification
<b>Isotype</b>	IgG
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SMARCA4

## Image



Western blot analysis of extracts of various cell lines, using SMARCA4 antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.