



# Recombinant Rabbit Elongation factor 1-delta (EEF1D)

<b>Product Code</b>	CSB-EP007431RB
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P53787
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	TTNFLVHEK IWFDKFKYDD AERSFYERMN GPVPGPSRQE NGATVILRDI ARARENIQKS LAGSSGPGAS SGPGGDHSEL AVRIASLEVE NQNLRGVVQD LQRAVSKLEA RLSALEKSSP THRASAPQTQ HVSPMRQVEP PAKKAAAPAE DDEDDIDL F GSDEEEDKEA ARLREERLRQ YAEKKARKPA LVAKSSILLD VKPWDEETDM ARLEACVRSV QLDGLVWGAS KLVPVGYGIR KLQIQCVVED DKVGTDLLEE EITKFEHVQ SVDIAAFNKI
<b>Source</b>	E.coli
<b>Target Names</b>	EEF1D
<b>Protein Names</b>	Recommended name: Elongation factor 1-delta Short name= EF-1-delta
<b>Expression Region</b>	2-280
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit functions as guanine nucleotide exchange factor. It is reported that this subunit interacts with HIV-1 Tat, and thus it represses the translation of host-cell, but not HIV-1, mRNAs. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.