



Recombinant *Drosophila erecta* 40S ribosomal protein SA (sta)

| | |
|--------------------------|---|
| Product Code | CSB-EP020485DLK-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | B3P9J3 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila erecta</i> (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Sequence | SGGLDILSL KEDDITKMLV ATTHLGSENV NfqMEQYVYK RRADGVNILN LGKTWEKLQL AARAI AIDN PSDIFVISSR PIGQRAVLKF AKYTDTTPIA GRFTPGAFTN QIQPAFREPR LLVVTDPNTD HQPIMEASYV NIPVIAFTNT DSPLRYIDIA IPCNNKSAHS IGLMWWLLAR EVLRLRGTIS RSVEWPVVVD LFFYRDPEEA EKEEAAAKEL LPAPKIEEAV DHPVEETT NW ADEVAAETVG GVEDWNEDTV KTSWGSDGQF |
| Source | <i>E. coli</i> |
| Target Names | sta |
| Protein Names | Recommended name: 40S ribosomal protein SA Alternative name(s): Protein stubarista |
| Expression Region | 2-270 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |