



Recombinant Human IST1 homolog (IST1)

Product Code	CSB-MP012171HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53990
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLGSGFKAER LRVNLRVIN RLKLEKKT ELAQKARKEI ADYLAAGKDE RARIRVEHII REDYLVEAME ILELYCDLLL ARFGLIQSMK ELDSGLAESV STLIWAAPRL QSEVAELKIV ADQLCAKYSK EYGKLCRTNQ IGTVNDRLMH KLSVEAPPKI LVERYLIEIA KNYNVPYEPD SVVMAEAPPG VETDLIDVGF TDDVKKGGPG RGGSGGFTAP VGGPDGTVPM PMPMPMSAN TPFSYPLPKG PSDFNGLPMG TYQAFNIHP PQIPATPPSY ESVDDINADK NISSAQIVGP GPKPEASAKL PSRPADNYDN FVLPELPSVP DTLPTASAGA STSASEDIDF DDLSRRFEEL KKKT
Source	Mammalian cell
Target Names	IST1
Protein Names	Recommended name: IST1 homolog Short name= hIST1 Alternative name(s): Putative MAPK-activating protein PM28
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.