



# Recombinant *Drosophila willistoni* Methylthioribose-1-phosphate isomerase (GK22419)

<b>Product Code</b>	CSB-EP014799DMQ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B4NG41
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila willistoni</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSLQSIKYNR GSLEILDQLL LPVQSKYVAV RGVEDGWKVI NKMQVRGAPA IAIVGCLSLA VEIYPEEFES KKSLRQEIEG KLNLYLSARP TAVNMKIAAD ELLSLANDLT KAEDVSV PDM KQRFLNATEA MLQKDIADNR AIGANGAQAI LKR VATT LGA AGTSGNGSVR VLTHCNTGSL ATAGYGTALG VIRNLHEL GK LEHVYCTETR PYNQGARLTA YELVHEKFPA TLVLDSMVAA LLRAKNVA AV VVGADRVAAN GDTANKIGTY QIAVVAKHHG VPFYVAAPLT SIDLQIPSGD HIIIEERPDR EMTHVGEHRI AAPGINCWNP AFDVTPASLI TGIITERGVF EPQELKEAIT KLLGL
<b>Source</b>	E.coli
<b>Target Names</b>	GK22419
<b>Protein Names</b>	Recommended name: Methylthioribose-1-phosphate isomerase Short name= M1Pi Short name= MTR-1-P isomerase EC= 5.3.1.23 Alternative name(s): S-methyl-5-thioribose-1-phosphate isomerase Translation initiation factor eIF-2B sub
<b>Expression Region</b>	1-365
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.