



Recombinant Helicobacter pylori Ferrochelataase (hemH)

Product Code	CSB-EP008579HUZ-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B2UUI4
Product Type	Recombinant Protein
Immunogen Species	Helicobacter pylori (strain Shi470)
Purity	>85% (SDS-PAGE)
Sequence	MNLINEKLNN LENNATKSPK EAVILLNMGG PNSLYEYGVF LKNMFDDPFI LTIKNNFMRK MVGKMIVNSR IEKSKKIYEK LGGKSPLTPI TFALTERLNK LDPSRFYTYA MRYTPPYASM VLQDLALKEV ESLVFFSMYP QYSSTTTLSS FNDAFNALKA LETFRPKVRV IERFYASKKL NEIILNTILS TLNNRKSQDF VLIFSVHGLP KSVIDAGDTY QQECEHHVNL LKELVQQKNI SFKEVLLSYQ SKLGPMKWLE PSTEELIEKH RKSHIIIYPL AFTIDNSETL YELDMQYRLM AERLAVKEYL VCPCLNDSIE FAKFIIELVN NLKSE
Source	E.coli
Target Names	hemH
Protein Names	Recommended name: Ferrochelataase EC= 4.99.1.1 Alternative name(s): Heme synthase Protoheme ferro-lyase
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.