



Recombinant Mouse FAS-associated factor 1 (Faf1)

Product Code	CSB-EP007963MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P54731
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MASNMDREMI LADFQACTGI ENIDEAITLL EQNNWDLVAA INGVIPQENG ILQSDFGGET MPGPTFDPAS HPAPASTPSS SAFRPVMPSR QIVERQPRML DFRVEYRDRN VDVVLEDSCT VGEIKQILEN ELQIPVPKML LKGWKTGDVE DSTVLKSLHL PKNNSLYVLT PDLPPSSSSS HAGALQESLN QNFMLIITHR EVQREYNLNF SGSSTVQEVK RNVYDLTSIP VRHQLWEGWP ASATDDSMCL AESGLSYPCH RLTVGRRRTSP VQTREQSEEQ STDVHMVSDS DGDDFEDASE FGVDDGEVFG MASSTLRKSP MMPENAENEG DALLQFTAEF SSRYGDCHPV FFIGSLEAAF QEAFYVKARD RKLAIYLHH DESVLTNVFC SQMLCAESIV SYLSQNFITW AWDLTKDNR ARFLTMCNRH FGSVIAQTIR TQKTDQFPLF LIIMGKRSSN EVLNVIQGNT TVDELMMRLM AAMEIFSAQQ QEDIKDEDER EARENKREQ DEAYRLSLEA DRAKREAHER EMAEQFRLEQ IRKEQEEERE AIRLSLEQAL PPEPKEENAE PVSKLRIRTP SGEFLERRFL ASNKLQIVFD FVASKGFPWD EFKLLSTFPR RDVTQLDPNK SLLEVNLFPPQ ETLFLQAKE
Source	E.coli
Target Names	Faf1
Protein Names	Recommended name: FAS-associated factor 1
Expression Region	1-649
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	Interaction of Fas ligand (TNFSF6) with the FAS antigen (TNFRSF6) mediates programmed cell death, also called apoptosis, in a number of organ systems. This protein binds to FAS antigen and can initiate apoptosis or enhance apoptosis initiated through FAS antigen. Initiation of apoptosis by This protein requires a ubiquitin-like domain but not the FAS-binding domain.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.