



Recombinant *Drosophila erecta* Ubiquitin-conjugating enzyme E2 S (GG18165)

Product Code	CSB-BP025477DLK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3NWW9
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSSQYSNVEN LSPQTIQVM RELQEMETTP PEGIKVLIN SDVTDIQALI DGPAGTPYAA GVFRVKLTN KDFPQTTPKA YFLTKIFHPN VAANGEICVN TLKKDWKPD L GIKHILLTIK CLLIVPNPES ALNEEAGKML LERYDDYSQR ARMMTEIHAQ PAKCGAGAHG DDKDDDGPST KKHAGLDKKL QDKKKEKLLK EKKRMLKRL
Source	Baculovirus
Target Names	GG18165
Protein Names	Recommended name: Ubiquitin-conjugating enzyme E2 S EC= 6.3.2.19 Alternative name(s): Ubiquitin carrier protein S Ubiquitin-protein ligase S
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.