



# Recombinant Human Serine protease HTRA4 (HTRA4)

<b>Product Code</b>	CSB-MP010904HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P83105
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	EKLHTQPSC PAVCQPTRCP ALPTCALGTT PVFDLCRCCR VCPAAEREVC GGAQQGPCAP GLQCLQPLRP GFPSTCGCPT LGGAVCGSDR RTYPSMCALR AENRAARRLG KVPAPVPVQWG NCGDTGTRSA GPLRRNYNFI AAVVEKVAPS VVHVQLWGRL LHGSRLVPVY SGSGFVISED GLIITNAHVV RNQQWIEVVL QNGARYEAVV KDIDLKDLA VIKIESNAEL PVLMLGRSSD LRAGEFVVAL GSPFSLQNTA TAGIVSTKQR GGKELGMKDS DMDYVQIDAT INYGNSGGPL VNLGDVIGV NSLRVTDGIS FAIPSDRVRQ FLAEYHEHQM KGKAFSNKKY LGLQMLSLTV PLSEELKMHY PDFPDVSSGV YVCKVVEGTA AQSSGLRDHD VIVNINGKPI TTTTDDVKAL DSDSLMAVL RGKDNLLLT IPETIN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	HTRA4
<b>Protein Names</b>	Recommended name: Probable serine protease HTRA4 EC= 3.4.21.-
<b>Expression Region</b>	32-476
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.