



# Recombinant Mouse Ras-related protein Rab-8A (Rab8a)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP019222MO   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P55258   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Mus musculus (Mouse)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MAKTYDYLFK LLLIGDSGVG KTCVLF RFSE DAFNSTFIST IGIDFKIRTI<br>ELDGKRIKLQ IWDTAGQERF RTITTAYYRG AMGIMLVYDI TNEKSFDNIR<br>NWIRNIEEHA SADVEKMILG NKCDVNDKRQ VSKERGEKLA LDYGIKFMET<br>SAKANINVEN AFFTLARDIK AKMDKKLEGN SPQGSSHGVK ITVEQQKRTS<br>FFRC  |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | Rab8a  |
| <b>Protein Names</b>     | Recommended name: Ras-related protein Rab-8A Alternative name(s):<br>Oncogene c-mel  |
| <b>Expression Region</b> | 1-204  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full length protein  |
| <b>Target Details</b>    | This protein is a member of the RAS superfamily which are small GTP/GDP-binding proteins with an average size of 200 amino acids. The RAS-related proteins of the RAB/YPT family may play a role in the transport of proteins from the endoplasmic reticulum to the Golgi and the plasma membrane. This protein shares 97%, 96%, and 51% similarity with the dog RAB8, mouse MEL, and mouse YPT1 proteins, respectively and contains the 4 GTP/GDP-binding sites that are present in all the RAS proteins. The putative effector-binding site of this protein is similar to that of the RAB/YPT proteins. However, this protein contains a C-terminal CAAX motif that is characteristic of many RAS superfamily members but which is not found in YPT1 and the majority of RAB proteins. Although this gene was isolated as a transforming gene from a melanoma cell line, no linkage between MEL and malignant melanoma has been demonstrable. This oncogene is located 800 kb distal to MY09B on chromosome 19p13.1. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final   |



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.