



Recombinant *Drosophila erecta* Phosphotriesterase-related protein (GG16392)

Product Code	CSB-EP018966DLK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3P1R1
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSTVQTVLGT ITPNLLGRTL THEHVALDFE HFYRPPPPDF ECELKAKISM STLGYVRLYP YSSKENVRFY DGEALEAAKK DVLLYKKHGG GSIVENSSYG LKRNLFFIVE LAKSTGVHFI AGTGHYIHAM QDASHASLSV EQMSDLYSKD IITGLQVNGK MVKCGFIGEV ASVYPIHDFE KNAIKAAGEI QEVLGCGVSM HPRVTKAPF EIMRLYLEAG GRADKCVMSH LDRTIFDIDE LLEFAKLG CY IQYDLFGTEC SFYQLNTSVD MISDGQRIDN LIKLIKEGLV DKLLMSHDIH TKHRLTSYGG HGYHHIHTNI LPRMFDRGVT LEQVEQMTVT NPANWLA FDP
Source	<i>E. coli</i>
Target Names	GG16392
Protein Names	Recommended name: Phosphotriesterase-related protein EC= 3.1.-.- Alternative name(s): Parathion hydrolase-related protein
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.